

FIGHTING FUNGALS Results in Time to

Make a Difference

With a robust fungal menu, differential testing, faster results, and anti-fungal therapy monitoring. Eurofins Viracor can help you better support your patients continuum of care.

Invasive fungal diseases (IFDs) comprise a serious and often life-threatening group of infections caused by fungi that affect tissues and organs. These diseases are becoming increasingly common, especially in people with weakened immune systems due to complications of organ transplantation and other comorbidities such as HIV/AIDS and cancer. IFDs can be serious and life threatening, especially in those who are immunocompromised. Early diagnosis and appropriate treatment are essential in managing these infections and improving patient outcomes.



Fungal pathogens are a major threat to public health as they are becoming increasingly common and resistant to treatment with only four classes of antifungal medicines currently available, and few candidates in the clinical pipeline. Most fungal pathogens lack rapid and sensitive diagnostics and those that exist are not widely available or affordable globally.

The invasive forms of these fungal infections often affect severely ill patients and those with significant underlying immune system related conditions. Populations at greatest risk of invasive fungal infections include those with cancer, HIV/AIDS, organ transplants, chronic respiratory disease, and post-primary tuberculosis infection.

Emerging evidence indicates that the incidence and geographic range of fungal diseases are both expanding worldwide due to global warming and the increase of international travel and trade.¹

~ The World Health Organization



Testing to individualize your approach for IFD patients



Aspergillosis:

Aspergillus is a common fungus that is found in the environment and can cause a variety of infections in the lungs, sinuses, and other organs. Aspergillosis is the term used to describe the different infections caused by this fungus. The most common form of aspergillosis is allergic bronchopulmonary aspergillosis, which affects people with asthma or cystic fibrosis. Invasive aspergillosis is a more serious infection that can occur in people with weakened immune systems. Symptoms of invasive aspergillosis include fever, cough, chest pain, and shortness of breath. **Panel includes:** *A. fumigatus, A. terreus, and Pan-Aspergillus*



8900 - Aspergillus Real-time PCR Panel 🕨



Candidiasis:

This yeast is a normal part of the human body, but it can cause infections in people who are immunocompromised. Candida is a type of fungus that lives on the skin and in the mouth, throat, and genital area. However, it can sometimes overgrow and cause infections in the body. Candidiasis is the term used to describe these infections. The most common form of candidiasis is thrush, which is an infection of the mouth and throat. Invasive candidiasis is a more serious infection that can occur in people with weakened immune systems. Symptoms of invasive candidiasis include fever and chills, and the infection can spread to the blood and other organs. Swab sampling of skin for C. auris is available by Real-time PCR.

Panel includes: Candida albicans, C. glabrata, C. tropicalis, C. krusei, C. parapsilosis & C. auris



33306 - Candida Real-time PCR Panel





Coccidiomycosis: (aka: Valley Fever)

Coccidioides immitis can cause a lung infection called Valley Fever. Coccidiomycosis is caused by inhaling spores of the fungus Coccidioides immitis, commonly found in the soil of certain parts of the southwestern United States, Mexico, and Central and South America. Most people who inhale these spores do not get sick, but some can develop flu-like symptoms, and in rare cases, the infection can spread to other parts of the body and become life-threatening.



3400 - OMEGA Coccidioides Antibody EIA 🕨



Cryptococcosis:

Cryptococcus is typically found in the soil and in bird droppings. It can cause infections in people with weakened immune systems, such as those with HIV/AIDS. Cryptococcal meningitis is a serious infection that can occur when the fungus enters the brain and spinal cord. Symptoms of cryptococcal meningitis include headache, fever, and confusion.



1700 - Fungitell® B-D-Glucan Assay







Histoplasmosis:

Histoplasma capsulatum is a fungus found in soil contaminated with bird or bat droppings. Histoplasmosis is caused by inhaling spores of the fungus. Most people who inhale these spores do not get sick, yet some will develop flu-like symptoms. In rare cases, the infection can spread to other parts of the body and become life-threatening.



30881/30882/33222/33251 - Histoplasma Galactomannan EIA



Mucormycosis:

Mucormycosis, also known as zygomycosis, is a rare but serious fungal infection caused by fungi of the order Mucorales. It is commonly observed in patients with uncontrolled diabetes, those undergoing chemotherapy, and those with immunosuppression. The fungus can invade blood vessels, leading to tissue ischemia and necrosis. Common clinical manifestations include fever, headache, facial pain, and swelling. Diagnosis is typically made by histopathology and culture of affected tissues. Treatment involves aggressive surgical debridement of affected tissues, along with antifungal therapy such as amphotericin B.



3200 - Mucorales Real-time PCR 🕨



Nocardiosis:

Nocardia is a genus of aerobic, gram-positive bacteria that can cause a variety of infections, including pulmonary infections, brain abscesses, and cutaneous infections. In immunocompromised patients, disseminated nocardiosis can occur. Diagnosis can be challenging, as the clinical presentation can be non-specific, and the organism may not be readily visible on Gram stain. Culture and molecular techniques are often necessary for definitive diagnosis. Treatment typically involves a combination of sulfonamides and other antimicrobial agents, depending on the severity and location of the infection.



6800 - Nocardia Real-time PCR Panel 🕨



Pneumocystis pneumonia:

Pneumocystis pneumonia (PCP) is a serious infection caused by the fungus Pneumocystis jirovecii. This fungus is found in the lungs and can cause pneumonia. This fungus can cause severe pneumonia, especially in immunosuppressed individuals such as those with HIV/AIDS. Diagnosis is typically made by identifying the organism in respiratory secretions or by histopathology. Treatment involves a combination of sulfamethoxazole and trimethoprim. Pneumocystis jirovecii was formerly classified as a protozoan but is now considered a fungus.



2000 - Pneumocystis jiroveci Real-time PCR Panel 🕨



IFD Individual Assays from Eurofins Viracor

| Individual Assays | Order Code(s) | Turnaround Time | Assay Range Available Sample Days Pe Types Days Pe | | Days Performed | NY State |
|--|---|--|---|---|--|------------------------------|
| Aspergillus Galactomannan | 2009 - BAL | Same day (within 8-12 hours from receipt of specimen) | <0.500 - Negative <0.500- Positive | BAL, Bronch Wash, CSF, Serum | Mon - Sat | NY State Approved |
| Coccidiodes EIA | 3410 - Serum | 3 days from sample receipt | ays from sample Positive, Intermediate, Negative Serum Mon - Fri | | NY State Approved Testing performed by Euro- fins DPT | |
| Fungitell [(1,3)-β-D- glucan] | 1703 - CSF 1709 - BAL 1710 - Serum 1726 - Bronch Wash | Same day (within 8-12 hours from receipt of specimen) | le day in 8-12 hours receipt of simen) Semi-Quant: <60 pg/mL Negative, 60-79 pg/mL Intermediate ≥0 pg/mL Positive Assay range Varies by Specimen Type CSF, BAL, Serum, Bronch Wash Mon - Sat | | NY State Approved | |
| Fungitell w/reflex to titer [(1,3)-β-D-glucan] | 351703 - CSF 351709 - BAL 351710 - Serum 351726 - Bronch Wash | Same Day (within 8-12 hours from receipt of specimen) | Semi-Quant: <60 pg/mL Negative, 60-79 pg/mL Intermediate ≥0 pg/mL Positive Assay range Varies by Specimen Type | CSF, BAL, Serum, Bronch Wash | Mon - Sat | Urine - NY State Approved |
| Histoplasma Galac- tomannan | 30881 - Serum 30882 - Urine 33222 - BAL 33251 - CSF | Same day (within 8-12 hours from receipt of specimen) | Urine/Serum: < 0.3 ng/mL - Negative, ≥ 0.3 ng/mL - Positive BAL: <0.73 ng/mL are consid- ered positive CSF: < 0.2 ng/ml are considered negative, ≥ 0.2 but <0.8 ng/mL are considered positive but below the limit of quantitation (LOQ), ≥ 0.8 ng/ mL are considered positive | 0.3 ng/mL 3 ng/mL - mL are consid- nl are ative, ≥ 0.2 but considered ow the limit of O(), ≥ 0.8 ng/ red positive | | NY State Approved |
| Mucorales by Real- time PCR | 3205 - Tissue 3209 - BAL | Same day (within 8-12 hours from receipt of specimen) | Detected/Not Detected | BAL, Tissue | Mon - Fri | NY State Approved |
| Nocardia by Real- time PCR | 6809 - BAL | Same day (within 8-12 hours from receipt of specimen) | BAL: Detected/Not Detected | BAL | Mon - Sat | NY State Approved |
| Pneumocystis by Real-time PCR | 2009 - BAL 2016 - Whole Blood 2019 - Trach Asp 2029 - Bronch Wash 2048 - Trach Wash | 8-12 hours from receipt of specimen (Same day) | Whole blood 213 - 1x108 copies/mL BAL, Trach Asp, Bronch Wash, & Trach Wash 84 - 1x108 copies/mL | d 213 - 1x108 Asp, Bronch ach Wash copies/mL BAL, Whole Blood, Trach Asp, Bronch Wash, Trach Wash | | NY State Approved |
| C. auris (Swab) | 33330 | Target 24 hours - after specimen receipt | Detected/Not Detected The LOD for C. auris was determined to be 600 cfu/mL | Skin swab | Mon - Sat | |

FIGHTING FUNGALS Results in Time to Make a Difference

View Our Fungal Test Menu

IFD Fungal Panels from Eurofins Viracor

| Panels | Order Code(s) | Turnaround Time | Assay Range | Available Sample Types | Days Performed | NY State |
|--|--|---|---|--|----------------|----------------------|
| Aspergillus Real-time PCR Panel | 8909 - BAL 8926 - Bronch Wash 8909 - Fresh frozen tissue 8906 - FFPE tissue | BAL & Bronch Wash Same day (within 12 - 24 hours from receipt of speci- men) Tissue: 2-5 business days (from receipt of specimen) | Detected/Not Detected: Pan- Aspergillus, with speciation for A. fumigatus and A. terreus | BAL, Bronch Wash, Fresh frozen tissue, FFPE tissue | Mon - Sat | NY State Approved |
| Candida Real-time PCR Panel | 33306 - Whole Blood | Whole blood same day (Within 24 hours from receipt of specimen) | Detected/Not Detected: Candida albicans, C. glabrata, C. tropicalis, C. krusei, & C. parapsiloisis, C. auris | Whole blood | Mon - Sat | NY State Approved |
| Fungal Plus Real- time PCR Profile | PFL8066 - BAL | 12 - 24 hours from receipt of specimen (same day) | See individual Panels/As- says for ranges. Consists of Aspergillus Real-time PCR Panel, Mucorales Real-time PCR, and Nocardia Real-time PCR assays. | BAL | Mon - Fri | NY State Approved |
| Fungal Plus II Real- time PCR Profile | PFL8006 - BAL | 12 - 24 hours from receipt of specimen (same day) | See individual Panels/As- says for ranges. Consists of Aspergillus Real-time PCR Panel, Mucorales Real-time PCR, Nocardia Real-time PCR, Pneumocystis Real-time PCR assays. | BAL | Mon - Fri | NY State Approved |

Raising the general awareness of the significance and impact of IFD on human health, in both the hospital and the community, is critical to understand the scale of the problem and to raise interest to help fight these devastating diseases.¹



🔅 eurofins

 Firacative C. Invasive fungal disease in humans: are we aware of the real impact? Mem Inst Oswaldo Cruz. 2020 Oct 9;115:e200430. doi: 10.1590/0074-02760200430. PMID: 33053052; PMCID: PMC7546207. <u>https://pubmed.ncbi.nlm.nih.gov/33053052/</u>





When timing is critical, can you afford not to know?

RISKS OF TREATMENT DELAY

Treat IFDs Faster, Reduce Mortality

Managing fungal infections in immunocompromised individuals can be challenging due to a lack of specific symptoms and increased drug resistance. A rapid and accurate diagnosis is key in reducing the risk of mortality due to delayed treatment.^{1,2}

| Antifungal Treatment Delay In days | 0 | 1 | 2 | <u><</u> 3 | |
|--|-----|-----|------------------------|--------------------------|--|
| Mortality | 15% | 24% | 37% | 40% | |
| Blood Culture 50% sensitivity | | 1-5 | 1-5 days for result | | |
| Other lab's Fungitell [®] testing ≥ 20% failure rate | | | 2-4 days | for result | |
| Eurofins Viracor Fungitell [®] testing < 6% failure rate | | | 8-12 hour for Serun | rs (M-Sat.) n samples | |



Antifungal Therapy, Drug Level Monitoring and Azoles

Once a fungal infection has been diagnosed, treatment can begin. Antifungal medications, such as fluconazole, amphotericin B, or voriconazole, are often used to treat fungal infections. The choice of medication will depend on the type of fungus causing the infection, as well as the patient's overall health and medical history.

Azoles

Azole drugs are a class of antifungal drugs that work by blocking the growth of fungi. They are effective against a wide range of fungi, including Candida, Aspergillus, and Cryptococcus.

Azole drugs can have side effects, such as liver damage, kidney damage, and bone marrow suppression. It is important to talk to your doctor about the risks and benefits of azole drugs before starting treatment.

Antifungal therapy is an important component in the management of IFDs. Drug level monitoring is important to ensure that therapeutic levels of the drug are achieved and maintained. Azole drugs such as voriconazole and Posaconazole are commonly used for the treatment of IFDs, and therapeutic drug monitoring is recommended to ensure optimal dosing and to prevent toxicity. Additionally, these drugs can interact with other medications, so it is important to monitor for drug-drug interactions and adjust doses as necessary.

8

Anti-fungal TDM Offering(s) from Eurofins Viracor

| Anti-Fungal TDM | Order Code(s) | Turnaround Time | Assay Range | Available Sample Types | Days Performed | Platform | NY State |
|------------------------------|-------------------------------|--|---------------|---------------------------|----------------|----------|----------------------|
| lsavuconazole (CRESEMBA®) | 4901 - Plasma 4910 - Serum | 8 - 12 hours from receipt of speci- men (same day) | 0.1-10 mcg/mL | Plasma, Serum | Mon - Sat | LC-MS/MS | NY State Approved |
| Itraconazole | 2801 - Plasma 2810 - Serum | 8 - 12 hours from receipt of speci- men (same day) | 0.1-10 mcg/mL | Plasma, Serum | Mon - Sat | LC-MS/MS | NY State Approved |
| Posaconazole | 4210 - Serum | 8 - 12 hours from receipt of speci- men (same day) | 0.1-10 mcg/mL | Serum | Mon - Sat | LC-MS/MS | NY State Approved |
| Voriconazole | 3301 - Plasma 3310 - Serum | 8 - 12 hours from receipt of speci- men (same day) | 0.1-10 mcg/mL | Plasma, Serum | Mon - Sat | LC-MS/MS | NY State Approved |



Unfortunately these [fungal] infections are pretty non-specific... there's not a very characteristic sign of fungal infection like you would see, for example, with Mpox (monkeypox virus), and only when we do a lot of testing [do] we find out you have an invasive fungal infection.

Dr. Luis Ostrosky, Chief of Infectious Disease and Epidemiology, UT Health and Memorial Hermann





With a robust fungal menu, differential testing, faster result turnaround time and diagnostics to monitor drug levels, Eurofins Viracor can help you better support your patients' continuum of care — from the onset of symptoms to managing dosage.



* Most fungal tests results the same day (within 8-12 hours) of specimen receipt.

FUNGAL TESTING For Your Patient's Continuum of Care

View Our Fungal Test Menu

Get Fast Accurate Results

Find out more at <u>eurofins-viracor.com</u> or contact us at <u>info@eurofins-viracor.com</u> or call 800-305-5198

About Viracor

With over 30 years of diagnostic expertise in infectious disease, immunology and allergy testing for immunocompromised and critical patients, Eurofins Viracor is passionate about delivering accurate, timely and actionable results, never losing sight of the connection between the testing it performs and the patients it serves.

Eurofins Viracor is a subsidiary of Eurofins Scientific (EUFI.PA), a global leader in bio-analytical testing, and one of the world leaders in genomic services. For more information, please visit <u>eurofins.com</u> and <u>eurofins-viracor.com</u>

References

- Firacative C. Invasive fungal disease in humans: are we aware of the real impact? Mem Inst Oswaldo Cruz. 2020 Oct 9;115:e200430. doi: 10.1590/0074-02760200430. PMID: 33053052; PMCID: PMC7546207.
- Kohler JR, Hube B, Puccia R, Casadevall A, Perfect JR. Fungi that infect humans. Microbiol Spectr. 2017;5(3).
- 3 WHO fungal priority pathogens list to guide research, development and public health actions.
- https://www.healthline.com/health/fungal-infection/antifungal 4.
- Assoc of Cape Cod, Fungitell Assay Instructions for Use. 2011 (Feb) Garey et al, Clin Infect Dis. 2006; 43:25-31. Karageorgopoulos et al. Clin Inf Dis. 2011; 52:750-770 5 Luong ML et al. Clin Inf Dis. 2012; 52 (10):1218-1226. Nguyen ML, Wissel MC et al. Clin Inf Dis. 2012; 54:1240-1248. 6
- Andres Pascual, Thierry Calandra, Saskia Bolay, Thierry Buclin, Jacques Bille, and Oscar Marchetti, Voriconazole Therapeutic Drug Monitoring in Patients with Invasive 7. Mycoses Improves Efficacy and Safety Outcomes, Clinical Infectious Diseases, 2008:46(2);204 by permission of Infectious Diseases Society of America.
- 8. Ascioglu S, Rex JH, de Pauw B, et al. Defining opportunistic invasive fungal infections in immunocompromised patients with cancer and hematopoietic stem cell transplants: an international consensus. Clin Infect Dis. 2002;34:7-14.2.
- Walsh TJ, Anaissie EJ, Denning DW, et al. Treatment of aspergillosis: clinical practice guidelines of the Infectious Diseases Society of America. Clin Infect Dis.2008; 9 46:327-60.3.
- 10. Pfeiffer CD, Fine JP, Safdar N. Diagnosis of invasive asapergillosis using a galactomannan assay: a meta-analysis. Clin Infect Dis. 2006;42:1417–27
- Information derived from the OMEGA Cocci Ab EIA Test Kit package insert (Immuno-Mycologics, Inc.).
- 12. Pappagianis D, Zimmer BL. Serology of coccidioidomycosis. Clin Microbiol Rev. 1990 Jul;3(3):247-68.
- 13. Bronnimann DA, Galgiani JN, Coccidioidomycosis, Eur J Clin Microbiol Infect Dis, 1989 May:8(5):466-473.
- 14. Galgiani, JN. Coccidioides immitis. in Mandell, GL, Bennett, JE, Dolin, R, Eds., Principles and Practice of Infectious Diseases, Fifth Edition. Churchill Livingstone, New York, NY. 2000;2746-57
- 15. Cox, M., G. Pesek, D. Phan, and G. Woods. 2008. Utility of a Histoplasma capsulatum Enzyme Immunoassay for Diagnosis of Disseminated Histoplasmosis and Correlation With Disease Activity. Arch. Pathol. Lab Med 132:1512.
 Gomez, B. L., J. I. Figueroa, A. J. Hamilton, B. L. Ortiz, M. A. Robledo, A. Restrepo, and R. J. Hay. 1997. Development of a novel antigen detection test for histoplasmosis.
- J. Clin. Microbiol. 35:2618-2622.
- 17. Kauffman, C. A. 2008. Diagnosis of histoplasmosis in immunosuppressed patients. Curr. Opin. Infect. Dis. 21:421-425.
- 18. Kauffman, C. A. 2007. Histoplasmosis: a clinical and laboratory update. Clin. Microbiol. Rev. 20:115-132.
- 19. Wheat, J. 1994. Histoplasmosis: recognition and treatment. Clin. Infect. Dis. 19 Suppl 1:S19-S27.
- Yeo, S. F. and B. Wong. 2002. Current status of nonculture methods for diagnosis of invasive fungal infections. Clin. Microbiol. Rev. 15:465-484.
 Cáceres DH, Samayoa BE, Medina NG, Tobón AM, Guzmán BJ, Mercado D, Restrepo A, Chiller T, Arathoon EE, Gómez BL. 2018. Multicenter validation of commercial antigenuria reagents to diagnose progressive disseminated histoplasmosis in people living with HIV/AIDS in two Latin American countries. J Clin Microbiol 56:e01959-17. https://doi.org/10.1128/JCM.01959-17
- 22. Product Insert, Clarus Histoplasma GM Enzyme Immunoassay. IMMY, 2019-05-08.
- 23. CDC Web site. http://www.cdc.gov/fungal/mucormycosis/definition.html
- Vitale RG, de Hoog GS, Schwarz P, et al. Antifungal susceptibility and phylogeny of opportunistic members of the order Mucorales. J Clin Microbiol. 2012 Jan;50(1);66-75. 24 25. Kasai M., Harrington SM, Francesconi A, et al. Detection of a molecular biomarker for Zygomycetes by quantitative PCR assays of plasma, bronchoalveolar lavage and lung tissue in a rabbit model of experimental pulmonary zygomycosis. J Clin Microbiol. 2008 Nov;46(11);3690-702.
- 26. Dannoui E, Schwarz P, Slany M, et al. Molecular detection and identification of Zygomycetes species from paraffin-embedded tissues in a murine model of disseminated zygomycosis: a collaborative European Society of Clinical Microbiology and Infectious Disease (ESCMID) fungal infection study group (EFISG) evaluation. J Clin Microbiol. 2010 June;48(6):2043-6.
- 27. Bernal-Martinzez L, Buitrago MJ, Castelli MV, Rodriquez-Tudela JL and Cuenca-Estrella M. Development of a single tube multiplex PCR to detect the most clinically relevant Mucormycetes species. Clin Microbiol Infect. 2013 Jan;19(1):E1-7
- Hammond SP, Bialek R. Milner DA, Petschnigg EM, Baden LR, Marty FM. Molecular methods to improve diagnosis and identification of mucormycosis. J Clin Microbiol. 2011 June;49(6):2151-3
- 29. Hata DJ, Buckwalter SP, Pritt BS, Roberts GD, Wengenack NL (2008) Real-time PCR method for detection of zygomycetes .J Clin Microbiol. 2008 Jul;46(7):2353-8. 30. Data on file at Viracor Eurofins.
- 31. Couble, Andrée, et al. "Direct Detection of Nocardia spp. in Clinical Samples by a Rapid Molecular Method." Journal of Clinical Microbiology, American Society for Microbiology, Apr. 2005, www.ncbi.nlm.nih.gov/pmc/articles/PMC1081390/.
- 32. Hosseini-Moghaddam, SM, et al. "Nocardiosis after Solid Organ Transplantation." Nocardiosis after Solid Organ Transplantation Infectious Disease and Antimicrobial Agents, www.antimicrobe.org/t28.asp.
- 33. Nocardia species (Nocardiosis). (n.d.). Retrieved January 22, 2018, from http://www.antimicrobe.org/b117.asp
- 34. "Nocardiosis." Centers for Disease Control and Prevention, Centers for Disease Control and Prevention, 24 Feb. 2017, www.cdc.gov/nocardiosis/index.html.
 35. Steinbrink, Julie, et al. "Nocardia Infections: Comparison of Manifestations and Outcomes in Immunocompromised (IC) and Non-Immunocompromised (Non-IC) Patients Open Forum Infectious Diseases | Oxford Academic." OUP Academic, Oxford University Press, 4 Oct. 2017, https://academic.oup.com/ofid/article/4/suppl_1/S708/4295716/ Nocardia-Infections-Comparison-of-Manifestations?searchresult=1
- 36. Wilson, John W. "Nocardiosis: Updates and Clinical Overview." Mayo Clinic Proceedings, Mayo Foundation, Apr. 2012, www.ncbi.nlm.nih.gov/pmc/articles/PMC3498414/. 37. Astellas Pharma US, Inc., CRESEMBA® Package Insert.
- 38. von Mach MA1, Burhenne J, Weilemann LS. 2006. Accumulation of the solvent vehicle sulphobutylether beta cyclodextrin sodium in critically ill patients treated with intravenous voriconazole under renal replacement therapy. BMC Clin Pharmacol. 2006 6:6.
- 39. Pascual A1, Calandra T, Bolay S, Buclin T, Bille J, Marchetti O. 2008. Voriconazole therapeutic drug monitoring in patients with invasive mycoses improves efficacy and safety outcomes. Clin Infect Dis. 46(2):201-11.
- 40. Andes D, Pascual A, Marchetti O. Antifungal Therapeutic Drug Monitoring: Established and Emerging Indications. Antimicrob Agents Chemother. 2009 Jan; 53(1):24–34.
- 41. Goodwin ML, Drew RH. Antifungal serum concentration monitoring: an update. J Antimicrob Chemother. 2008 Jan; 61(1):17-25
- 42. Pappas PG1, Kauffman CA, Andes D, et al. Clinical practice guidelines for the management of candidiasis: 2009 update by the Infectious Diseases Society of America. Clin Infect Dis. 2009 Mar 1; 48(5):503-35.
- 43. Walsh TJ1, Anaissie EJ, Denning DW, et al. Treatment of aspergillosis: clinical practice guidelines of the Infectious Diseases Society of America. Clin Infect Dis. 2008 Feb 1; 46(3):326-60
- 44 Ashbee HR, Barnes RA, Johnson EM, et al. Therapeutic drug monitoring (TDM) of antifungal agents: guidelines from the British Society for Medical Mycology. J Antimicrob Chemother 2014; 69: 1162-76
- 45. Pascual A, Calandra T, Bolay S, et al. Voriconazole therapeutic drug monitoring in patients with invasive mycoses improves efficacy and safety outcomes. Clin Infect Dis 2008; 46: 201-11.
- 46. Tan K, Brayshaw N, Tomaszewski K, et al. Investigation of the potential relationships between plasma voriconazole concentrations and visual adverse events or liver function test abnormalities. J Clin Pharmacol 2006; 46: 235-43.







Eurofins Viracor 18000 W 99th St., Ste. #10 • Lenexa, KS 66219 Ph 800.305.5198 • Fax 816.347.0143 • Eurofins-Viracor.com

MM1452B CC23033B